ATLAS Diabetes Self-test

Ref. 70004001

Self-test that measures the presence of glucose and/or ketones in urine.

INTENDED USE:

The Atlas Diabetic Check kit contains dip-and-read test strips and is intended for use to check for glucose and ketones levels in urine specimens as an aid in diagnosing diabetes. The test provides results by the visual comparison with a colour chart printed on the pouch. Each colour denotes a range of concentrations of glucose and ketones.

INTRODUCTION:

Diabetes is a common disease caused by high sugar (glucose) levels in blood. This condition is called Hyperglycaemia, and is caused by reduced insulin secretion or action. Insulin is a hormone that is secreted from the pancreas. The major role of this hormone is to regulate sugar levels in the blood.

There are two types of diabetes: I and II. Type I (known as Insulin-Dependent type) is characterised by absolute deficiency of insulin production. Type II(Non-Insulin-Dependent type) is characterised by relative deficiency of insulin, which means that insulin is produced in lower amounts or it does not act properly. As the level of glucose increases in the blood above the kidney's capacity it starts to appear in the urine. Therefore, measuring the presence of glucose in the urine can be used as tool to check for diabetes.

People with diabetes are more susceptible to developing specific complications such as retinopathy (a pathological disorder of the retina) blindness, kidney failure, nerve damage and circulatory problems, which may lead to heart disease and stroke. Therefore detection of glucose and seeking early medical assistance as applicable could help to prevent the development of these symptoms.

In diagnosis of diabetes, doctors always advise to check for ketones to further confirm the diabetic status. Ketones are groups of organic compounds that result from the metabolism of fatty acids. Non-diabetic people use glucose as an energy source, but in diabetic people, glucose is not passed to the cells, so cells will use fatty acids as the energy source instead.

The less glucose absorbed by cells (as in cases when the diabetes is more advanced) the more fatty acids are consumed and accordingly more ketone is produced. This condition is known as diabetic ketoacidosis. As the level of ketone increases in the blood, it will start to appear in urine.

Contents:

- Diabetes test
- Drying agent (do not use)
- Manual
- Color chart

PRECAUTIONS:

1. Please read all the information in this leaflet before performing the test.
2. Do not use the test after the expiration date.
3. If the package is not completely sealed do not use the test.
4. Do not open the test foil pouch until it has reached room temperature and you are ready to start the test.
5. The test should be performed in a well-lit area.
6. Use the test device immediately after opening it.
7. Do not touch the test area. This could affect results and may also be harmful.
8. Use a disposable sample container to be discarded after performing the test.
9. The pouch contains a silica gel pack to absorb humidity. Do not open the pack. Throw it away with the remainder of the test.
10. Do not freeze.
11. At the end of the test, wrap everything you have used in a plastic bag and throw in the bin. Do not forget to wash your hands properly.
12. The remaining sample should be discarded and flushed in the toilet.
13. Keep out of the reach of children.
14. For in vitro diagnostic and self-testing use. Not to be taken internally.

STORAGE:

- Store at room temperature between 15°-30° (59°F-86°F).
- Do not store the strips in the refrigerator or freezer.
- Do not expose strips to moisture, heat and light before use.
- Use the strip immediately after removing it from the pouch.


SPECIMEN COLLECTION AND PREPARATION:

Collect fresh urine sample in a clean and dry disposable container. Make sure there are no traces of detergent on the container. Test the urine as soon as possible after collection.

PROCEDURE:

This procedure **MUST BE FOLLOWED EXACTLY** to achieve reliable test results.

1. Check that the product is within the expiration date shown on the kit pack.
2. Prepare the urine specimen.
3. Remove the strip from the pouch. Familiarise yourself with the position of the reaction area of both ketone and glucose. The white reaction area is for ketones and the green area is for glucose. Also, familiarise yourself with the colour chart on the pouch.

4. Dip the test strip in the urine until the reaction areas are completely immersed for no more than 1 second.
5. Remove the dipstick from the urine and tap the strip on the rim of the cup to remove excess urine and place it horizontally with the reaction areas facing up.
6. Leave the strip for 30-60 seconds for the reaction to take place.
7. Read the results by comparing the colors of the reaction on the strip to those of the chart. While comparing, keep the strip in a horizontal position to avoid possible mixing of colors between the two reaction areas on the strip.
8. Identify the best color match on the color chart and the correspondent concentration range. A change in color that appears only along the edges of the reaction areas indicates that the reaction did not take place properly so we recommend redoing the test with another strip. Results read after 60 seconds are not valid.

MANUAL SELF-TEST ATLAS Diabetes Self-test

EN

quicktest.eu

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EN

Glucose	Ketones	Result
Negative	Negative	Not diabetic
Trace	Negative	Repeat the test. If this result continues to appear, consult your doctor. You may be developing diabetes.
+ / ++ / +++ / ++++	Negative	You could be diabetic. Seek medical assistance.
Negative	Trace	Not diabetic. The presence of ketones is not significant.
Trace	Trace	Repeat the test. If this result continues to appear, consult your doctor. You may be developing diabetes.
+ / ++ / +++ / ++++	Trace	You could be diabetic. You may be developing a condition called diabetic ketoacidosis. Seek urgent medical assistance.
Negative	+	Not Diabetic. The presence of Ketones may be associated with the presence of Other conditions such as diarrhea, vomiting or weight loss.
Trace	+	Repeat the test. If this result continues to appear, consult your physician. You may be developing Diabetes.
+ / ++ / +++ / ++++	+	You could be diabetic. You may be developing a condition called diabetic ketoacidosis. Seek urgent medical assistance.

Read results

EN

Glucose	Ketones	Result
Negative	++	Not diabetic. The presence of Ketones may be associated with the presence of other conditions such as diarrhea, vomiting or weight loss.
Trace	++	Repeat the test. If this result continues to appear, consult your doctor. You may be developing diabetes.
+ / ++ / +++ / ++++	++	You could be diabetic. Most probably you are developing a condition called diabetic ketoacidosis. Seek urgent medical assistance.
Negative	+++ / ++++	Not diabetic. The presence of ketones may be associated with the presence of Other conditions such as diarrhoea, vomiting or weight loss.
Trace	+++ / ++++	Repeat the test. If this result continues to appear, consult your doctor. You may be developing diabetes, although-the presence of high ketones is not likely with low glucose levels.
+ / ++ / +++ / ++++	+++ / ++++	You could be diabetic. Most probably you are developing a condition called diabetic ketoacidosis. Seek urgent medical assistance.

Questions and answers

EN

When is the best time to do the test?

This test can be done anytime during the day. It is good to do the test once before eating and once 2 hours after a meal.

If I dipped the strip for more than 1 second, is the result affected?

The reaction takes place immediately after dipping the strip. Leaving the strip for 2-3 seconds will not affect the result, but if kept for more than 5 seconds, results cannot be considered accurate and so you need to repeat the test.

I left the strip for 10 minutes and then I read the results, shall I consider the results true?

Results read after more than 60 seconds (1 minute) must not be considered correct so you need to repeat the test.

What are the signs of diabetes?

The most frequent signs of diabetes are: being very thirsty, urinating often, feeling very hungry or tired, losing weight, having sores that heal slowly, having dry itchy skin, having blurry eyes, losing the feeling in your feet or having tingling in your feet.

What could I do to control glucose level?

Doing the following will certainly help controlling glucose level: follow a diet that is advised by your doctor or dietician, be active and do some exercises for around 30 minutes each day and stop smoking.

ATLAS MEDICAL
William James House, Cowley Road, Cambridge, CB4 0WX, UK
Tel: +44 (0) 1223 858 910, Fax: +44 (0) 1223 858 524

	Read instructions before use		For in-vitro diagnostics		Do not reuse
	Store at +15° - +30°C		Batch number		Expiry date
	Manufacturer	PPI137A01 Rev F (21.04.2015)			

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