Urinary Tract Infection Self-test

Ref. 70021001
Zelftest die witte bloedcellen, nitrieten en/of bloed in de urine detecteert.

INTENDED USE

The Atlas Urinary Tract Infection Test provides dip-and-read test strips that are intended for use to check for Leukocytes, Nitrate and Blood in urine specimens as an aid in the diagnosis of urinary tract infection. The test provides results by visual comparison with a color chart printed on the pack.

Important! Read instructions carefully before use.

Kit components:

- 2 urinary tract infection tests
- 1 drying agent (not used)
- 1 manual
- 1 color chart

Precautions

1. Please read all the information in this leaflet before performing the test.
2. Do not use the test after the expiration date.
3. If the package is not completely sealed do not use the test.
4. Do not open the test foil pouch until it has reached room temperature and you are ready to start the test.
5. The test should be performed in a well-lit area.
6. Use the test device immediately after opening it.
7. Do not touch the test area. This could affect results and may also impose personal hazards.
8. Use a disposable sample container to be discarded after performing the test.
9. The pouch contains a Silica Gel pack to absorb humidity. Do not open the pack. Throw it away with the remaining of the test.
10. Do not freeze.
11. At the end of the test, wrap every thing you have used in a plastic bag and throw away in the bin. Do not forget to wash your hands properly.

Result interpretation

Leukocytes	Nitrite	Blood	Result
Positive	Positive	Positive	You have a urinary tract infection. There is presence of bacteria. The presence of blood is associated with severe infections. Seek medical assistance immediately.
		Negative	You have a urinary tract infection. There is presence of bacteria. Seek medical assistance immediately.
Negative	Positive	Positive	You have a urinary tract infection. You might have bacteria that cannot produce nitrite. Presence of blood is associated with severe infections. Seek medical assistance immediately.
		Negative	You have a urinary tract infection. You might have bacteria that cannot produce nitrite. Seek medical assistance immediately.

EN

MANUAL SELF-TEST URINARY TRACT INFECTION

EN

quicktest.eu

LIMITATIONS OF THE TEST:

Substances that cause abnormal urine color, such as some drugs may affect the color development on the strip. The color development on the reagent pad may be masked, or a color reaction may be produced on the pad that could be interpreted visually as a false positive. It is therefore recommended that in case of doubt, the test should be repeated after stopping the medication.

Leukocytes: The test result may not always be consistent with the leukocyte cell number by the microscopic examination. Positive results may be found with high humidity and high temperature condition, and failure of the bottle security. Positive results may occasionally be found in random specimens from females due to contamination of the specimens by the vaginal discharge.

Nitrite: Any degree of uniform pink color development should be considered positive; however, pink spots or pink edges should not be interpreted as a positive result. Color development is not proportional to the number of bacteria present. The urine test detects only nitrite-producing bacteria. Occasionally bacteria will be present that do not produce nitrite. Therefore, a negative result cannot rule out the presence of urinary tract infection.

Blood: A false positive can sometimes occur when bacteria are present in the urine. Ascorbic acid or protein may reduce the reactivity of the blood test. Strong oxidizing substances, such as hypochlorites, may produce a false positive result. Urine from menstruating females often, but not always, yields positive results.

HOW TO DETERMINE POSITIVE OR NEGATIVE VALUES:

Any color other than the color indicating negative result is considered positive. Refer to the results chart.

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FR

12. The remaining sample should be discarded and flushed in the toilet.
13. Keep out of the reach of children.
14. For in vitro diagnostic and self-testing use. Not to be taken internally.

STORAGE

- Store at room temperature between 15°-30° (59°F-86°F).
- Do not store the strips in the refrigerator or freezer.
- Do not expose strips to moisture, heat and light before use.
- Use the strip immediately after removing it from the pouch.

SPECIMEN COLLECTION AND PREPARATION

Collect a clean catch urine sample by washing the genital area and collect a midstream sample of urine in a clean and dry disposable container. The container has to be devoid of any detergent traces. This method of collecting urine helps prevent bacteria around the genital area from getting into the sample and confusing the test results. Test the urine as soon as possible after collection. A first morning urine specimen is recommended for nitrite test.

Result interpretation

Leukocytes	Nitrite	Blood	Result
Positive	Positive	Positive	This result is highly unlikely to occur and we recommend redoing the test. If you have the same result again, this means that you most probably have a urinary tract infection since there is presence of bacteria. Presence of blood is associated with severe infections. Seek medical assistance immediately.
		Negative	This result is highly unlikely to occur and we recommend redoing the test. If you have the same result, this means that you have a urinary tract infection as there is presence of bacteria. Seek medical assistance immediately.
Negative	Positive	Positive	No urinary tract infection. The presence of blood may be associated with other conditions such as kidney stones, excessive exercise, prostate diseases (in males) and some kidney or blood diseases. Seek medical assistance.
		Negative	No urinary tract infection.

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FR

PROCEDURE: This procedure MUST BE FOLLOWED EXACTLY to achieve reliable test results.

1. Check that the product is within the expiration date shown on the pack.
2. Prepare the urine specimen.
3. Remove the strip from the pouch. Familiarize yourself with the position of the reaction area of Leukocytes, RBCs and Nitrite. Pale pink reaction area is for Leukocytes, white is for nitrite and dark yellow area is for blood. Also, familiarize yourself with the color chart on the pack.
4. Dip the test strip in the urine until the reaction areas are completely immersed for no more than 1 second.
5. Remove the dipstick from the urine and tap the strip on the rim of the cup to remove excess urine and place it horizontally with the reaction areas facing up.
6. Leave the strip for 30-60 seconds for the reaction to take place.
7. Read the results by comparing the colors of the reaction on the strip with those of the chart. While comparing, keep the strip in a horizontal position to avoid possible mix of colors between the reaction areas on the strip.
8. Identify the best match color on the color chart and the correspondent concentration range. A change in color that appears only along the edges of the reaction areas indicates that the reaction did not take place properly so we recommend redoing the test with another strip. Results read after 60 seconds are not valid.

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